

Date _____

Customer code C- _____

Invoice address _____

Contact person _____

Your reference _____

Other delivery address _____

☐ Invoice address = delivery address☐ Quote

DECON E-DRIVE



Please add order form for the standard chair.

Mounting on the Netti wheelchair is included – the standard wheels will remain with the manufacturer.

☒ = Standard equipment

| Item no. | | Description | Price in € |
|----------------------------|-------------------------------------|--|------------------|
| D70001 | <input checked="" type="checkbox"/> | Decon E-Drive with standard configuration as stated below | 7.077,- |
| D71016 | <input type="checkbox"/> | Mounting kit for Netti chairs (required) | 208,- |
| D70013 | <input type="checkbox"/> | Delivery with standard Netti mail wheels for wheel combination 24" + 7" incl. quick release and bushers (a constant change between Decon and Netti wheels is only possible as special adaptation) | 74,- |
| D800 | <input type="checkbox"/> | Delivery with standard Netti mail wheels for other wheel combinations than 24" + 7" incl. quick release, bushers and other necessary parts (a constant change between Decon and Netti wheels is only possible as special adaptation) | Price on request |
| Battery and charger | | | |
| Standard | <input checked="" type="checkbox"/> | NiMH battery with a max range of approx. 18 km | 0,- |
| D70011 | <input type="checkbox"/> | Li-ion battery with a max range of approx. 32 km | 851,- |
| Standard | <input checked="" type="checkbox"/> | Battery charger 100-240V | 0,- |
| D70021 | <input type="checkbox"/> | Power converter for car 12V | 152,- |
| User weight | | | |
| Standard | <input checked="" type="checkbox"/> | Max user weight 125 kg | 0,- |

| | | | |
|--------|--------------------------|---|-------|
| D70031 | <input type="checkbox"/> | Heavy duty version - max user weight 150 kg by 20", 22", 24" or 26" wheels Rims and spokes are made from a stronger material | 730,- |
|--------|--------------------------|---|-------|

| Wheel size and accessories | | | |
|----------------------------|-------------------------------------|--|-------|
| D71015 | <input type="checkbox"/> | Anti-tippers (required - except if combined with a universal device holder) | 231,- |
| D71000 | <input type="checkbox"/> | 16" Flexel fine tread wheels (max. capacity 100 kg) | 326,- |
| D71001 | <input type="checkbox"/> | 20" fine tread pneumatic wheels with aluminium push rims | 326,- |
| D71002 | <input type="checkbox"/> | 22" fine tread pneumatic wheels with aluminium push rims | 326,- |
| Standard | <input checked="" type="checkbox"/> | 24" x 1 3/8 fine tread pneumatic wheels with aluminium push rims | 0,- |
| D71003 | <input type="checkbox"/> | 26" x 1" fine tread pneumatic wheels with aluminium push rims (only for the Netti III family) | 357,- |
| D71004 | <input type="checkbox"/> | PU wheels x1 (available for 20", 22", 24" and 26" x 1") | 92,- |
| D71014 | <input type="checkbox"/> | PU wheels wide x1 3/8 (available for 22" and 24" x 1 3/8") | 362,- |
| D71005 | <input type="checkbox"/> | Spoke guards, transparent (available for 20"-26" wheels) | 231,- |



[A]



[B]

| Driving unit | | | |
|---------------------|-------------------------------------|---|-------|
| Standard/ D71006 | <input checked="" type="checkbox"/> | Joystick for user, swing-away arm [A] <input checked="" type="checkbox"/> Right <input type="checkbox"/> Left | 0,- |
| D71007R/ D71007L | <input type="checkbox"/> | Joystick mounted for use by caregiver <input checked="" type="checkbox"/> Right <input type="checkbox"/> Left | 0,- |
| D71008R | <input type="checkbox"/> | Attendant control instead of joystick [B] <input checked="" type="checkbox"/> Right | 0,- |
| D71009R/ D71009L | <input type="checkbox"/> | Second joysticks for caregiver [2xA] <input checked="" type="checkbox"/> Right <input type="checkbox"/> Left | 866,- |
| D71010R | <input type="checkbox"/> | Attendant control in addition to the joystick [A+B] <input checked="" type="checkbox"/> Right | 830,- |



| Mounting | | |
|----------|-------------------------------------|--|
| Standard | <input checked="" type="checkbox"/> | Behind the back of the wheelchair, including bag for mounting |
| SPECIAL | <input type="checkbox"/> | For chairs with dynamic back unit a special holder for the battery under the frame must be requested |

Please note:

- Netti 4U chairs adjusted in seat depth 500 mm require frame extension.
- The mounting of wide front wheels (min. 45 mm) is required.
- For chairs with seat width extension, please contact customer service.

Specifications

| | |
|--|---|
| Motor | Yamaha motors: 24V 120W |
| Battery types | NiMH battery, 24V x 6.7Ah - range 18 km Li-ion battery, 25V x 11.8Ah - range 32 km |
| Wheel | 7 kg |
| Total weight | 16 kg (including NiMH battery) 17.5 kg (including Li-ion battery) |
| Total chair width | Netti 4U Family: Seat width approx. + 230 mm Netti III Family without frame extension: Seat width + 245 mm Netti III Family with frame extension*: Seat width +256 mm Netti III Family with 16" Decon E-Drive: Seat width approx. + 260 mm E-Drive increases the total width with approx. 40 mm *Please note that frame extension is included in some options like Extra-long seat depth, dynamic back solutions etc. Please contact customer service for questions. |
| Wheel dimensions | 16", 20", 22", 24", 26" |
| Air pressure | 6 Bar, 600 kPa |
| Speed | 0-6 km/t |
| Maximum range with fully charged battery | Approx. 18 km with NiMH battery/ Approx. 32 km with Li-ion battery in accordance with ISO 7176-4 |

| | |
|------------------------------|--|
| Fuse | 30 A (flat) |
| Max user weight | Up to 125 kg (for 16'': up to 100 kg) |
| Max user weight "Heavy-duty" | Up to 150 kg Rims and spokes are made from a stronger material |
| Max gradient | 6° |
| CE mark | Complies with EU-Directive 93/42/EEC |
| Joystick | Contactless sensing by single-chip/ Std. fitted on right side swing-away arm. |
| Force alteration | By optional spring (for individual optimization) |
| USB charger | Auxiliary USB port (5V, 500 mA max.) |
| Display device | LCD (Vertical aligned) |
| Parameter setting | Onboard (joystick) / via PC (Smart Tune) |
| Battery charger | 100 to 240V AC, 50/60Hz (standard) A power converter for car is available as accessory. |
| Rated power | 37V. 26A (when charging) |
| Charging time | Normal charging time approx. 2-3 hours |
| Possible combinations | Please send us a request for a combination with other options of the selected wheelchair as some combinations might not be possible. |