

Netti Precision seat cushion



Date _____

Invoice address _____

Customer code C- _____

Contact person _____


Reference _____

Other delivery address _____

☐ Invoice address = delivery address

☐ Quote

Netti Precision 1		Netti Precision 2	
			
	• When user is at risk of skin and pressure injuries. • Pressure distribution and relief. • Protects bony areas. • Support of the pelvic ensures a good seating position. • Great immersion minimises risk of sliding forward.		• For users with sensitive skin, high risk of skin damage or even skin injuries. • Enhanced pressure relief of tail bone and pelvic floor. • Pressure distribution and relieve. • Support of the pelvic ensures a good seating position. • Great immersion minimises risk of sliding forward

Different choices of covers	
 EASY CARE+ Durable, breathable, elastic, wipeable and hygienic. PES (Polyester) + Polyurethane.	 STELLA+ Breathable, elastic and durable. Enhances the microclimate, ventilates and reduces build-up of heat. This strengthens the skin and the resistance of the underlying tissue against internal and external influences, such as pressure and shear. The material is polyester and elastane.

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post.klepp@Meyragroup.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info.ry@Meyragroup.com | Phone +45 87887300

Netti Precision seat cushion



Netti Precision 1					
Cushion with Easy Care+ cover	Cushion with Stella+ cover	Width	Depth	Height	
☐ 234303	☐ 232303	350 mm	400 mm	65 mm	
☐ 234305	☐ 232305	350 mm	450 mm	65 mm	
☐ 234307	☐ 232307	350 mm	500 mm	65 mm	
☐ 234318	☐ 232318	380 mm	400 mm	65 mm	
☐ 234320	☐ 232320	380 mm	450 mm	65 mm	
☐ 234321	☐ 232321	380 mm	500 mm	65 mm	
☐ 234326	☐ 232326	400 mm	400 mm	65 mm	
☐ 234328	☐ 232328	400 mm	450 mm	65 mm	
☐ 234330	☐ 232330	400 mm	500 mm	65 mm	
☐ 234332	☐ 232332	400 mm	550 mm	65 mm	
☐ 234343	☐ 232343	430 mm	400 mm	65 mm	
☐ 234345	☐ 232345	430 mm	450 mm	65 mm	
☐ 234347	☐ 232347	430 mm	500 mm	65 mm	
☐ 234349	☐ 232349	430 mm	550 mm	65 mm	
☐ 234353	☐ 232353	450 mm	400 mm	65 mm	
☐ 234355	☐ 232355	450 mm	450 mm	65 mm	
☐ 234357	☐ 232357	450 mm	500 mm	65 mm	
☐ 234359	☐ 232359	450 mm	550 mm	65 mm	
☐ 244301	☐ 242301	500 mm	450 mm	65 mm	
☐ 244302	☐ 242302	500 mm	500 mm	65 mm	
☐ 244304	☐ 242304	500 mm	550 mm	65 mm	
☐ 244323	☐ 242323	550 mm	500 mm	65 mm	
☐ 244325	☐ 242325	550 mm	550 mm	65 mm	
☐ 244343	☐ 242343	600 mm	500 mm	65 mm	
☐ 244345	☐ 242345	600 mm	550 mm	65 mm	

☐ Other sizes available on request:

Precision 1 Width _____ Depth _____ Cover: ☐ Stella+ ☐ Easy Care+

Delivery time for Netti Precision seat cushions is 2-3 weeks.

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post.klepp@Meyragroup.com | Phone +47 51 78 62 20
 DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info.ry@Meyragroup.com | Phone +45 87887300

Netti Precision seat cushion



Netti Precision 2					
Cushion with Easy Care+ cover	Cushion with Stella+ cover	Width	Depth	Height	
☐ 264303	☐ 262303	350 mm	400 mm	75 mm	
☐ 264305	☐ 262305	350 mm	450 mm	75 mm	
☐ 264307	☐ 262307	350 mm	500 mm	75 mm	
☐ 264318	☐ 262318	380 mm	400 mm	75 mm	
☐ 264320	☐ 262320	380 mm	450 mm	75 mm	
☐ 264321	☐ 262321	380 mm	500 mm	75 mm	
☐ 264326	☐ 262326	400 mm	400 mm	75 mm	
☐ 264328	☐ 262328	400 mm	450 mm	75 mm	
☐ 264330	☐ 262330	400 mm	500 mm	75 mm	
☐ 264332	☐ 262332	400 mm	550 mm	75 mm	
☐ 264343	☐ 262343	430 mm	400 mm	75 mm	
☐ 264345	☐ 262345	430 mm	450 mm	75 mm	
☐ 264347	☐ 262347	430 mm	500 mm	75 mm	
☐ 264349	☐ 262349	430 mm	550 mm	75 mm	
☐ 264353	☐ 262353	450 mm	400 mm	75 mm	
☐ 264355	☐ 262355	450 mm	450 mm	75 mm	
☐ 264357	☐ 262357	450 mm	500 mm	75 mm	
☐ 264359	☐ 262359	450 mm	550 mm	75 mm	
☐ 274301	☐ 272301	500 mm	450 mm	75 mm	
☐ 274302	☐ 272302	500 mm	500 mm	75 mm	
☐ 274304	☐ 272304	500 mm	550 mm	75 mm	
☐ 274323	☐ 272323	550 mm	500 mm	75 mm	
☐ 274325	☐ 272325	550 mm	550 mm	75 mm	
☐ 274343	☐ 272343	600 mm	500 mm	75 mm	
☐ 274345	☐ 272345	600 mm	550 mm	75 mm	

☐ Other sizes available on request:

Precision 2 Width _____ Depth _____ Cover: ☐ Stella+ ☐ Easy Care+

Delivery time for Netti Precision seat cushions is 2-3 weeks.

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post.klepp@Meyragroup.com | Phone +47 51 78 62 20
 DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info.ry@Meyragroup.com | Phone +45 87887300