

# HANDLE WITH CARE

## **Lock/Unlock**

Use the lever(s) located on the wheelchair base, front or back, to lock/unlock wheels. Use caution as this will allow chair to roll freely. Do not force or push the chair when in lock mode.

## **Lifting the chair**

Only carry from the base. Wheelchair is bottom heavy.

## **Keep the chair standing upright**

Do not place chair on it's side. Sensitive components can be damaged. Always keep chair standing upright.

## **Fuse to be switch off**

Ensure that main fuse is switched off. Ask owner of chair for fuse's location.

## **Bring wheelchair to cabin door.**

For questions, always consult the owner of the chair:

Contact name: \_\_\_\_\_

Phone number: \_\_\_\_\_

Wheelchair model: \_\_\_\_\_

Measurements:      Weight: \_\_\_\_\_ Height: \_\_\_\_\_ Width: \_\_\_\_\_

Battery type: \_\_\_\_\_

Please refer to user manual for additional safety tips and transportation guidelines.