

Product Combination Agreement

Alu Rehab A.S

and

Symmetric Designs

Bedriftsvegen 23
4353 Klepp Stasjon
Norway

125 Knot Place
Salt Spring Island, BC
Canada, V8K 2M4

As manufacturers of

Wheelchairs:

Netti 4U CE, Netti 4U CE Pluss, Netti 4U Base
Netti 4U CED and CED Dynamic, Netti S,
Netti Dynamic S, Netti III, Netti III HD, Netti III el.
Netti Dynamic III HD, Netti Dynamic AdaptPro
Netti AdaptPro

Head supports and fixing:

Savant P, S and A head positioning
Axion Rotary Interface 6025 and 9025
Netti Adaptor NA 150 and NA 250,
Savant SP, for fixing head supports

Netti A and E head supports and fixing kit 89558

hereby declare that the above-mentioned products can be combined under the following conditions:

- 1) The Netti fixing kit 89558 is constructed for fixing Netti Head supports type A and type E to the Axion Rotary Interface 6025 and Axion Rotary Interface 9025.
- 2) The Netti Adaptors NA-150 and NA 250 are constructed for being fixed the Axion Rotary Interfaces 6025 and 9025 and fitting into the profile on Netti Head support vertical bar with latch lever art. No 41402.
- 3) The Savant head support can be fixed to the Axion Rotary Interface and the Netti Adaptor which enables them to fit in the profile of the Netti vertical pole with latch and thereby being fixed to all listed Netti wheelchairs.
- 4) Netti head supports A and E fixed to Axion Rotary Interface which again can be fixed to the Netti Adaptors allows for the head supports to be fixed to Netti wheelchairs allowing the head supports to move freely in a rotary movement.
- 5) It is the responsibility of the person fitting the wheelchair head support to the user to secure that the adjustments of the head support are done in a way that serves the user most beneficial.
- 6) The user instruction for each product must be read carefully.
The user instruction contains important information and warnings to make the products to be safe to use for the user.
- 7) The user instruction for the use of the Savant head support and the Axion Rotary Interface is established by Symmetric Designs.
- 8) The user instruction for the use of Netti wheelchairs are established by Alu Rehab.
- 9) Wheelchairs fitted with Savant head support and / or Axion Rotary Interface have not been crash tested and are not approved for being used as a seat in a vehicle. If a crash tested Netti

wheelchair is to be used as a seat in a vehicle the Savant head support and head supports with Axion Rotary Interface has to be replaced by a standard Netti Head support during transport.

- 10) This agreement between *Alu Rehab AS* and Symmetric Designs (the Parties) states that the combination of the *Savant head support, the Axion Rotary Interface and the Netti Adaptor* and *Alu Rehab AS* Head support bracket for Netti wheelchairs in combination (Products) meets the Medical Device Regulation MDR 2017/745 (EU) concerning product combination.
- 11) The Products are both separately CE-marked devices and registered separately at a notified body.
- 12) The Parties are responsible for possible limitations in use, instructions for installation and use, risk analysis, and required documentation of the medical device(s) listed above. *Alu Rehab AS* pledges to inform *Symmetric Designs* about any design changes to its wheelchairs which may affect the safety and or function of the Products in combination. *Symmetric Designs* pledges to inform *Alu Rehab A.S* about any design changes to its head support devices which may affect the safety and or function of the Products in combination.
- 13) Each Party has the product liability for its medical device. The Parties jointly have the product liability for the Products in combination.
- 14) The Parties renounce responsibility for a combination if not assembled according to each Parties installation and use instructions.

Should any incident occur with the combination of the products described in this agreement, a copy of this agreement should be added to the report to the Authorities.

Provided the above conditions are met, the manufacturers hereby declare that the combination of the products also meet the requirements according to the standards.

This agreement is valid from today and is subject to three months' notice for termination.

This document has been drawn up in two copies of which each Party has taken one.

If regulatory requirements should change or if constructional changes are made to the products, each of the parties shall evaluate whether or not this affects the validity of this contract and inform the other party without any delay of such circumstances.

This agreement shall be valid for an indefinite period of time as from this execution. However it can be terminated by either of the parties by 30 days given notice in writing.

Klepp Stasjon / Salt Spring Island

Dato February 2022

Alu Rehab AS



Björn Carlzon
CEO

Dato February 2022

Symmetric Designs

Samuel Hannah
Sam Hannah
CEO