


## Monopad description (not available for Netti 4U CE or CE Plus)

The monopad will be mounted in the tubes where the extension pieces for the upper leg support normally are mounted. Only removable with tools.

If mounted on existing chair, please state serial number: \_\_\_\_\_

If mounted on new chair, please attach order form.



The angle of the monopad can be regulated with gas spring from 90-45°. If a fixed angle of the knee position is wanted, please state the wished angle (adjustable with tools): \_\_\_\_\_

Leg length measured from seat plate to footboard (standard is 400-530 mm): \_\_\_\_\_

### Calf support specifications

Calf support width (max width = seat width of the chair): \_\_\_\_\_

Calf support height (max height = 220 mm): \_\_\_\_\_

Calf support thickness (standard is 25 mm): \_\_\_\_\_

Calf support cover if other than Skai (3D or Easy Care available): \_\_\_\_\_

**Footboard aluminium specifications**

The footboard is foldable upwards. The angle of the footboard can be adjusted 10° downwards with tools. At the back there is a cut out to prevent possible conflict with the front wheels.

Footboard width front (max width = seat width of the chair): \_\_\_\_\_

Footboard width back (standard is app. seat width minus 160 mm): \_\_\_\_\_

Footboard depth (max depth = 350 mm): \_\_\_\_\_

Need for slits for the mounting of ankle huggers or other foot fixations: \_\_\_\_\_